Die intrauterine Insemination (IUI),

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Intrauterine insemination (IUI),

in-vitro fertilisation with embryo transfer (IVF/ET) and related procedures in assisted reproduction basically bear the same risks of cross-contamination for the patient and a resulting pregnancy as would be the case in natural conception.

However, in reproductive-medical intervention involved doctors adopt special responsibility for the health of the treated woman as well as of the child which could result from the treatment.

AWMF Guidelines 014/040 (2002)

Technical data

Length:	230 mm
External diameter:	2,60 mm
Internal diameter:	2,20 mm
Internal canula diameter:	1,10 mm
Catheter volume	0,35 ml

Delivery in sterile packaging. Application only if packaging is undamaged up to the printed expiry date.

**For single use only!
Do not sterilise!**

Insemination catheter, 10 per package

Order No.: MED 1000003
PZN: 0661374
EAN: 4250201200051

Insemination catheter, 1 per package

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SP 00498.000 RRN 21030

Inseminations-Katheter Insemination-Catheter

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Insemination-Catheter

The medesign insemination catheter is a special medical product specifically and exclusively developed for intrauterine insemination treatment (IUI). The medesign catheter is particularly suitable for intrauterine insemination because its flexibility enables easy insertion through the cervix while simultaneously providing excellent stability for further insertion.

Delivered as a disposable, ready-for-use device in sterile packaging.

Intrauterine insemination treatment is indicated, if the ascension of a sufficient number of mobile and morphologically normal sperm into the upper female genital tract is not guaranteed, e.g.

- Oligo-Asthenozoospermia
- post-surgical stenoses of the cervical canal (e.g. after conisation) with consecutive dysmucorrhea
- negative post-coital test acc. to Sims-Huhner

Principles of treatment

Aided by the medesign insemination catheter, sperm is injected directly into the vagina. In order to gain a sufficient number of mobile sperm and to remove substances from the ejaculate which could lead to cramping of the uterus or infections, the ejaculate requires previous professional preparation. Suitable procedures for sperm preparation are Swim-Up or density gradient centrifugation. They also provide an effective physical germ reduction. Principally, the best possible sterility must be observed for IUI to prevent crosscontamination in the uterus: The use of sterile materials, wearing hygienic clothing, sterile gloves and mouth protection is essential.

Application of the medesign insemination catheter

Reclined on a gynaecological treatment chair, the cervix is set with a (self-adjusting) speculum in such a way, that the opening of the cervix is clearly visible. The sterile medesign insemination catheter is fastened to an insulin syringe, whose plunger has been retracted to 0.5 ml.

Check secure and firm seating of the connection! Draw up 0.4 ml cultured medium of rinsed, washed sperm via the catheter. The catheter should be filled with liquid.

Use sterile swab on cervix and cervical neck. Grasp the medesign insemination catheter with corn pliers approx. 2-3 cm from catheter tip. Strictly sterile rinsing of soft, flexible medesign insemination catheter through the cervix into the uterus, preferably without touching the vaginal walls.

Advance to approx. 4-5 cm up to Fundus Uteri. At stop retract the catheter by approx. 5mm. Push the wings attached the handle to the height of the tubes. (The openings discharging the liquid are located at the same height as the wings and are aimed directly at the tubes.)

Use gentle, even pressure on the plunger of the syringe and slowly inject the sperm medium into the Cavum uteri.

The air at the end of the syringe serves to completely empty the catheter.

With slight rotating movement remove the catheter from the uterus slowly after approx. 60 - 90 seconds.

In order to keep the sperm in the uterus, the inseminated woman should remain in a slightly reclined position for a further 20 minutes (pelvis tilted by approx. 30°) on the examination chair. The use of a portio cap is neither required nor sensible.

Application by medical staff only! Remarks on intrauterine insemination

Approx. 30 minutes prior to the intrauterine insemination, the patient rectally receives a spasmolyticum as suppository to suppress any occurring uterus spasms caused by the irritating sperm concentrate.

Approx. 20 hours after insemination, LH control is carried out on the urine and/or an ultrasound scan of the follicle. Is the follicle still intact and LH values have risen, a second IUI will be conducted.

The probability of pregnancy per IUI cycle is approx. 10 - 15 %. This however is only attained if the IUI is carried out shortly before ovulation!